

RICH-MAR CORPORATION

THERAMINI 2

Preventive Maintenance Procedure

Serial No. _____

Date _____

- (1) Power up the unit while depressing the Stop/Clear button, look for the Rich-Mar logo verify software version. 1.03

Quadpolar Function Tests

Using a test load 500 ohm 5 watts to a oscilloscope.

Set up the scope to read 50 uS/division and 20 volts/division.

Set up TMT 1 as Quadpolar, 80-150 Hz, Vector - shallow/normal, Continuous, 10:00.

Start TMT 1. Increase intensity to the maximum.

- (1) Verify the peak to peak voltage of channel 1 on your oscilloscope is approximately 70 volts.
- (2) Verify the peak to peak voltage of channel 2 is approximately 70 volts.

Monophasic Function Tests

Connect Channel 1 to the oscilloscope. Set the scope to read 50 volts/division.

Set up TMT 2 as Monophasic, 100Hz, 50 - 50 uS, Continuous, 10:00.

Start TMT 2. Increase intensity to the maximum on both channels.

- (1) Verify that the peak voltage on Channel 1 is 100 volts.
- (2) Verify that the peak voltage on channel 2 is 100 volts.
- (3) Change the scope to 1 mS/division. See that the twin pulses repeat every 10 divisions.

BiPhasic Function Tests

Connect channel 1 to the oscilloscope. Set the scope to read 50 Volts/division and 50 uS/division.

Set up TMT 3 as Biphasic, 100 Hz, 200 - 50 uS, surge -10 "ON" and 10 "OFF", 10:00.

Start TMT 3. Increase the intensity to the maximum.

- (1) Verify that the peak to peak voltage on channel 1 and channel 2 is 200.
- (2) Verify that the unit is "ON" for 10 seconds and "OFF" for 10 seconds on both channels.

BiPolar Function Tests

Connect channel 1 to the oscilloscope. Set the scope to read 10 V/division and 50 uS/division.

Set up TMT 4 as bipolar, 100 Hz, med/slow, surge - 10 "ON" and 10 "OFF", 10:00.

Start TMT 4. Increase intensity to the maximum.

- (1) View the wave form on your oscilloscope. The wave form should vector between 60 and 80 V peak to peak at a rate of 10 cycles in 30 seconds on channel 1 and channel 2.
- (2) Verify that the output is "ON" for 10 seconds and "OFF" for 10 seconds.

Russian Function Tests

Connect channel 1 to the oscilloscope. Set the scope to read 20 volts/division and 100 uS/division

Set up TMT 5 as Russian. 50 Hz, surge - 15 seconds "ON" and 50 seconds "OFF", 10:00.

Start TMT 5. Increase intensity to maximum.

- (1) View the wave form on your oscilloscope. Verify that it repeats itself every 4 divisions (both channels)
- (2) Verify that the peak to peak voltage is 6 divisions on both channels.
- (3) Verify that the output is "ON" for 15 seconds and "OFF" for 50 seconds.

Micro Function Tests

Connect channel 1 to the oscilloscope. Set the scope to read 500mV/division and 500 uS/division.

Set up TMT 6 as micro, 500 Hz, Continuous, 10:00.

Start TMT 6. Increase the intensity to the maximum.

- (1) Verify that the peak voltage is approximately 2 divisions on both channels.
- (2) Verify that the peaks are 4 divisions apart on both channels.
- (3) Verify that the output alternates from positive to negative at approximately 2.7 second intervals.

Safety and Appearance Tests

- (1) Check the units operation at 108 and 132 line volts AC.
- (2) Check the line leakage using a Simpson 229-2 line leakage meter.
Forward_____ Reflected _____
- (3) Verify that the patient lead cord test functions correctly.
- (4) Recheck the overall appearance of the unit.
- (5) Check encoder pot for proper function when tured screen should move in same direction .
- (6) Check panel for proper touch control on all settings.

All checks completed by _____